### FILLING LINE FOR APTAR CSP VIALS




# FILLING LINE FOR APTAR CSP VIALS

- PROJECT THAT NEEDS TO BE ADJUSTED AFTER PRE-SERIES
- FOR NORTHERN AMERICA ONLY
- TAILOR-MADE SYSTEM, ADJUSTABLE TO FIT VARIOUS REQUIREMENTS
- SPEED FROM 15 TO 60 VIALS/MIN, DEPENDING ON THE OPTIONS SELECTED, AND THE VIALS AND TABLETS TO PACK
- CHANCE PARTS TO MANAGE A WIDE RANGE OF VIALS
- A BETTER PRODUCTIVITY
- OUR EXPERTISE AT YOUR SERVICE
- OUR SUPPORT DURING IMPLEMENTATION

Elizabeth Europe, expert tooling partner for pharmaceutical industry. The European Production site of the Elizabeth Companies is located in France. Elizabeth Europe has a turnover of 12 million euros, has a staff of 100 people, and has 5 activities : Compression tooling (1<sup>st</sup> in France, 650 customers in Europe), Blister tooling (1<sup>st</sup> in France, 300 customers in Europe), Special Machineries, Press parts (1<sup>st</sup> in France, 50 % of French market) and Eliza-Press (rotary press).



## C2 FILLING LINE FOR APTAR CSP VIALS



## • A. VIAL FILLING SYSTEM

Purpose:

It allows faster automatic opening and closing of the vials before and after filling. This system bypasses all handling.



Features:

### - All size of vials

- Specific but adjustment to vial diameter 17,5mm/height 40,8 mm
- Hand supply of the empty vials into the unit A
- Maximum machine speed is 60 vials/min
- But usual speed is 15 to 20 vials/min due to hand supply of the vial
- Data management through password functions
- Emptiness
- Proper closing checking
- Remote access to automation
- Change parts specific to vial diameter
- CFR 21 part 11 optional

## FILLING LINE FOR APTAR CSP VIALS

### • **B.** AUTOMATED VIAL SUPPLY UNIT

#### Purpose:

Unit dedicated to automatic supply of the empty vials into unit A.

### Features:

- Maximum speed of unit B is 60 vials per min

- The final speed depends on vial size, tablet size/surface state/ number, cotton and opening/closing features


- All size of CSP vials

- But vial diameter 17,5 mm and height 40,8 mm may require specific adjustment (not included)

- Change parts specific to vial diameter

- Aluminium, stainless steel and food compliant plastic

- Vials loaded manually into a bowl feeder located close to the ground, easily accessible to operator



### • C. INTEGRATED COUNTING UNIT

Purpose:

Unit dedicated to a high-speed automatic counting of tablets/capsules and filling into the vials.

#### Features:

- 12 tracks counting unit, one funnel

- The counting sensors are capable of automatically adjusting against light reduction due to dust. So, the sensors are dust friendly and ensures count accuracy ever for very dusty tablets. - If the dust blocks the sensors, a warning signal appears on


the monitor.

- Handle coated and uncoated tablets and hard/soft gel capsules

Ex: For tablet diameter 5 mm, speed 216.000 tablets/hour Ex: For tablet diameter 10 mm, speed 117.000 tablets/hour

- Ex: For capsule size 0, speed 63.000 capspules/hour
- Ex: For capsule size 2, speed 90.000 capspules/hour
- Adjustable vibrator for tablet flow

- In option (not included), sorting device by air purging and vacuum suction to remove broken tablets



Vial	Capsule	Count	Speed in Vial/min (+/- 5%)
M-3003-200 64 ml	Size 0	30	50
M-3003-200 64 ml	Size 1	45	60
M-3003-200 64 ml	Size 2	60	46
M-3003-477 110 ml	Size 0	45	48
M-3003-477 110 ml	Size 1	60	54
M-3003-477 110 ml	Size 2	90	34
M-3003-451 166 ml	Size 0	90	30
M-3003-451 166 ml	Size 1	120	38
M-3003-451 166 ml	Size 2	150	18

## FILLING LINE FOR APTAR CSP VIALS

	<b>Price Target</b> min <sub>-</sub> max		
A- Semi automated vial filling unit			
A.1 Equipment			
Mechanical base, 1 format, FAT, IQOQ doc, packaging	\$60 000	-	\$70 000
A.2 Options & customization			
Additional format (tube diameter)	\$1 500	-	\$2 500
Additional closing/opening device	\$500	-	\$1 000
Special parts, special components	quotation	unde	er request
Link with another counting unit	quotation under request		
A.3 Services recommended			
fixed cost travel From Europe to the USA (2 technicians)	\$5 000	-	\$5 800
commissioning (2 days, 2 technicians)	\$4 000	-	\$4 600
SAT (1 day, 2 technicians)	\$2 000	-	\$2 300
Training (1 day, 2 technicians)	\$2 000	-	\$2 300
Assistance to manufacturing start (1 day x 2 technicians)	\$2 000	-	\$2 300
B- Automated vial supply unit			
B.1 Equipment			
Mechanical base, 1 format, FAT, IQOQ doc, packaging	\$80 000	-	\$90 000
B.2 Options & customization			
additional format (tube diameter)	\$500	-	\$1 000
Special parts, special components	quotation	unde	er request
B.3 Services recommended			
commissioning (2 days, 2 technicians)	\$2 000	-	\$2 300
SAT (1 day, 2 technicians)	\$2 000	-	\$2 300
Training (1 day, 2 technicians)	\$1 000	-	\$1 200
Assistance to manufacturing start (1 day x 2 technicians)	\$1 000	-	\$1 200
C- Integrated Counting unit			
C.1 Equipment			
Mechanical base, 1 format, FAT, IQOQ doc, packaging	\$115 000	-	\$125 000
C.2 Options & customization			
additional format (funnel)	\$500	-	\$1 000
Special parts, special components	quotation	unde	er request
C.3 Services recommended			
commissioning (2 days, 2 technicians)	\$500	-	\$600
SAT (1 day, 2 technicians)	\$500	-	\$600
Training (1 day, 2 technicians)	\$500	-	\$600
Assistance to manufacturing start (1 day x 2 technicians)	\$500	-	\$600

Due to continued updating of our equipment, the texts, illustrations and technical specifications mentioned on these pages are illustrative only and are subject to change without notice.

**ELIZABETH EUROPE** ZI Les Gailletrous 1 Rue Mickaël Faraday 41260 LA CHAUSSÉE SAINT VICTOR - FRANCE Ph. : +33 (0)2 54 90 21 20 Fax : +33 (0)2 54 74 75 26 E-mail : service@elizeurope.com www.elizeurope.com

